



LMAN1 Antibody, Biotin conjugated

Product Code	CSB-PA012991LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49257
Immunogen	Recombinant Human Protein ERGIC-53 protein (154-329AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein ERGIC-53 (ER-Golgi intermediate compartment 53 kDa protein) (Gp58) (Intracellular mannose-specific lectin MR60) (Lectin mannose-binding 1), LMAN1, ERGIC53 F5F8D
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LMAN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.