



LRP2 Antibody, Biotin conjugated

Product Code	CSB-PA013096LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P98164
Immunogen	Recombinant Human Low-density lipoprotein receptor-related protein 2 protein (1211-1369AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Low-density lipoprotein receptor-related protein 2 (LRP-2) (Glycoprotein 330) (gp330) (Megalin), LRP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LRP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.