



MAP1A Antibody, Biotin conjugated

Product Code	CSB-PA013400LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78559
Immunogen	Recombinant Human Microtubule-associated protein 1A protein (1979-2168AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Microtubule-associated protein 1A (MAP-1A) (Proliferation-related protein p80) [Cleaved into: MAP1A heavy chain; MAP1 light chain LC2], MAP1A, MAP1L
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.