



MEIS3 Antibody, Biotin conjugated

Product Code	CSB-PA013689LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99687
Immunogen	Recombinant Human Homeobox protein Meis3 protein (178-264AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein Meis3 (Meis1-related protein 2), MEIS3, MRG2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MEIS3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.