



METAP1 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA013715LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P53582 |
| Immunogen | Recombinant Human Methionine aminopeptidase 1 protein (54-125AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Methionine aminopeptidase 1 (MAP 1) (MetAP 1) (EC 3.4.11.18) (Peptidase M 1), METAP1, KIAA0094 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | METAP1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |