



METAP1 Antibody, Biotin conjugated

Product Code	CSB-PA013715LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53582
Immunogen	Recombinant Human Methionine aminopeptidase 1 protein (54-125AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Methionine aminopeptidase 1 (MAP 1) (MetAP 1) (EC 3.4.11.18) (Peptidase M 1), METAP1, KIAA0094
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	METAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.