



NEK4 Antibody, HRP conjugated

Product Code	CSB-PA015702LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51957
Immunogen	Recombinant Human Serine/threonine-protein kinase Nek4 protein (516-661AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase Nek4 (EC 2.7.11.1) (Never in mitosis A-related kinase 4) (NimA-related protein kinase 4) (Serine/threonine-protein kinase 2) (Serine/threonine-protein kinase NRK2), NEK4, STK2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NEK4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.