




NID1 Antibody

Product Code	CSB-PA015802LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14543
Immunogen	Recombinant Human Nidogen-1 protein (792-942AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:200-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nidogen-1 (NID-1) (Entactin), NID1, NID
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NID1

Image



IHC image of CSB-PA015802LA01HU diluted at 1:200 and staining in paraffin-embedded human melanoma performed on a Leica Bond™ system.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.