



NPHP4 Antibody, Biotin conjugated

Product Code	CSB-PA015987LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75161
Immunogen	Recombinant Human Nephrocystin-4 protein (397-543AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nephrocystin-4 (Nephroretinin), NPHP4, KIAA0673
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NPHP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.