



NRDC Antibody, HRP conjugated

Product Code	CSB-PA016075LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43847
Immunogen	Recombinant Human Nardilysin protein (32-113AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nardilysin (EC 3.4.24.61) (N-arginine dibasic convertase) (NRD convertase) (NRD-C) (Nardilysin convertase), NRDC, NRD1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NRDC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.