



PAM Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA017417LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P19021 |
| Immunogen | Recombinant Human Peptidyl-glycine alpha-amidating monooxygenase protein (338-497AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Peptidyl-glycine alpha-amidating monooxygenase (PAM) [Includes: Peptidylglycine alpha-hydroxylating monooxygenase (PHM) (EC 1.14.17.3); Peptidyl-alpha-hydroxyglycine alpha-amidating lyase (EC 4.3.2.5) (Peptidylamidoglycolate lyase) (PAL)], PAM |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | PAM |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |