



PAM Antibody, Biotin conjugated

Product Code	CSB-PA017417LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19021
Immunogen	Recombinant Human Peptidyl-glycine alpha-amidating monooxygenase protein (338-497AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidyl-glycine alpha-amidating monooxygenase (PAM) [Includes: Peptidylglycine alpha-hydroxylating monooxygenase (PHM) (EC 1.14.17.3); Peptidyl-alpha-hydroxyglycine alpha-amidating lyase (EC 4.3.2.5) (Peptidylamidoglycolate lyase) (PAL)], PAM
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PAM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.