




PDLIM4 Antibody

Product Code	CSB-PA017735LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50479
Immunogen	Recombinant Human PDZ and LIM domain protein 4 protein (81-243AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IP; Recommended dilution: WB:1:500-1:5000, IP:1:200-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PDZ and LIM domain protein 4 (LIM protein RIL) (Reversion-induced LIM protein), PDLIM4, RIL
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PDLIM4

Image



Western Blot

Positive WB detected in: K562 whole cell lysate


All lanes: PDLIM4 antibody at 6.1µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 36, 26 kDa

Observed band size: 36 kDa



Immunoprecipitating PDLIM4 in K562 whole cell lysate

Lane 1: Rabbit control IgG instead of CSB-PA017735LA01HU in K562 whole cell lysate. For western blotting, a HRP-conjugated Protein G antibody was used as the secondary antibody (1/2000)

Lane 2: CSB-PA017735LA01HU (6 μ g) + K562 whole cell lysate (500 μ g)

Lane 3: K562 whole cell lysate (20 μ g)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.