



PDLIM4 Antibody, HRP conjugated

Product Code	CSB-PA017735LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50479
Immunogen	Recombinant Human PDZ and LIM domain protein 4 protein (81-243AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PDZ and LIM domain protein 4 (LIM protein RIL) (Reversion-induced LIM protein), PDLIM4, RIL
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PDLIM4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.