



PDLIM4 Antibody, FITC conjugated

Product Code	CSB-PA017735LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50479
Immunogen	Recombinant Human PDZ and LIM domain protein 4 protein (81-243AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PDZ and LIM domain protein 4 (LIM protein RIL) (Reversion-induced LIM protein), PDLIM4, RIL
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PDLIM4