



PHOX2A Antibody, FITC conjugated

Product Code	CSB-PA017940LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14813
Immunogen	Recombinant Human Paired mesoderm homeobox protein 2A protein (150-264AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Paired mesoderm homeobox protein 2A (ARIX1 homeodomain protein) (Aristaless homeobox protein homolog) (Paired-like homeobox 2A), PHOX2A, ARIX PMX2A
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PHOX2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.