



PRCP Antibody, HRP conjugated

Product Code	CSB-PA018636LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42785
Immunogen	Recombinant Human Lysosomal Pro-X carboxypeptidase protein (347-443AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysosomal Pro-X carboxypeptidase (EC 3.4.16.2) (Angiotensinase C) (Lysosomal carboxypeptidase C) (Proline carboxypeptidase) (Prolylcarboxypeptidase) (PRCP), PRCP, PCP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PRCP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.