



PROK1 Antibody, Biotin conjugated

Product Code	CSB-PA018746LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P58294
Immunogen	Recombinant Human Prokineticin-1 protein (20-105AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prokineticin-1 (Endocrine-gland-derived vascular endothelial growth factor) (EG-VEGF) (Mambakine), PROK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PROK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.