



PSAT1 Antibody, Biotin conjugated

Product Code	CSB-PA018838LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y617
Immunogen	Recombinant Human Phosphoserine aminotransferase protein (104-206AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphoserine aminotransferase (EC 2.6.1.52) (Phosphohydroxythreonine aminotransferase) (PSAT), PSAT1, PSA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PSAT1