



RBFOX2 Antibody, FITC conjugated

Product Code	CSB-PA019439LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43251
Immunogen	Recombinant Human RNA binding protein fox-1 homolog 2 protein (8-110AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA binding protein fox-1 homolog 2 (Fox-1 homolog B) (Hexaribonucleotide-binding protein 2) (RNA-binding motif protein 9) (RNA-binding protein 9) (Repressor of tamoxifen transcriptional activity), RBFOX2, FOX2 HRNBP2 RBM9 RTA
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RBFOX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.