



RFX2 Antibody, Biotin conjugated

Product Code	CSB-PA019611LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48378
Immunogen	Recombinant Human DNA-binding protein RFX2 protein (1-130AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-binding protein RFX2 (Regulatory factor X 2), RFX2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RFX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.