



RGS6 Antibody, Biotin conjugated

Product Code	CSB-PA019658LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49758
Immunogen	Recombinant Human Regulator of G-protein signaling 6 protein (177-262AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regulator of G-protein signaling 6 (RGS6) (S914), RGS6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RGS6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.