



RNF39 Antibody, FITC conjugated

Product Code	CSB-PA019888LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2S5
Immunogen	Recombinant Human RING finger protein 39 protein (200-321AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RING finger protein 39 (Protein HZFw), RNF39, HZFW
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RNF39
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.