



SCAMP1 Antibody, FITC conjugated

Product Code	CSB-PA020750LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15126
Immunogen	Recombinant Human Secretory carrier-associated membrane protein 1 protein (2-134AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Secretory carrier-associated membrane protein 1 (Secretory carrier membrane protein 1), SCAMP1, SCAMP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SCAMP1