



SERPINB4 Antibody, HRP conjugated

Product Code	CSB-PA021072LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48594
Immunogen	Recombinant Human Serpin B4 protein (128-208AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serpin B4 (Leupin) (Peptidase inhibitor 11) (PI-11) (Squamous cell carcinoma antigen 2) (SCCA-2), SERPINB4, PI11 SCCA2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SERPINB4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.