



SLC7A2 Antibody, Biotin conjugated

Product Code	CSB-PA021714LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52569
Immunogen	Recombinant Human Cationic amino acid transporter 2 protein (604-658AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cationic amino acid transporter 2 (CAT-2) (CAT2) (Low affinity cationic amino acid transporter 2) (Solute carrier family 7 member 2), SLC7A2, ATRC2 CAT2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC7A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.