



TAL1 Antibody, HRP conjugated

Product Code	CSB-PA023110LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17542
Immunogen	Recombinant Human T-cell acute lymphocytic leukemia protein 1 protein (1-114AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell acute lymphocytic leukemia protein 1 (TAL-1) (Class A basic helix-loop-helix protein 17) (bHLHa17) (Stem cell protein) (T-cell leukemia/lymphoma protein 5), TAL1, BHLHA17 SCL TCL5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TAL1