



TAL1 Antibody, Biotin conjugated

Product Code	CSB-PA023110LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17542
Immunogen	Recombinant Human T-cell acute lymphocytic leukemia protein 1 protein (1-114AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell acute lymphocytic leukemia protein 1 (TAL-1) (Class A basic helix-loop-helix protein 17) (bHLHa17) (Stem cell protein) (T-cell leukemia/lymphoma protein 5), TAL1, BHLHA17 SCL TCL5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TAL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.