



TOX Antibody, Biotin conjugated

Product Code	CSB-PA024073LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94900
Immunogen	Recombinant Human Thymocyte selection-associated high mobility group box protein TOX protein (141-288AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thymocyte selection-associated high mobility group box protein TOX (Thymus high mobility group box protein TOX), TOX, KIAA0808
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TOX