



TRHR Antibody, Biotin conjugated

Product Code	CSB-PA024441LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P34981
Immunogen	Recombinant Human Thyrotropin-releasing hormone receptor protein (320-398AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thyrotropin-releasing hormone receptor (TRH-R) (Thyroliberin receptor), TRHR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TRHR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.