



ZFP36 Antibody, Biotin conjugated

Product Code	CSB-PA026448LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26651
Immunogen	Recombinant Human mRNA decay activator protein ZFP36 protein (16-103AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	mRNA decay activator protein ZFP36 (G0/G1 switch regulatory protein 24) (Growth factor-inducible nuclear protein NUP475) (Tristetraprolin) (Zinc finger protein 36) (Zfp-36), ZFP36, G0S24 NUP475 RNF162A TIS11A TTP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ZFP36
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.