



MT3 Antibody, HRP conjugated

Product Code	CSB-PA04639B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25713
Immunogen	Recombinant Human Metallothionein-3 protein (1-68AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metallothionein-3 (MT-3) (GIFB) (GIF) (Growth inhibitory factor) (Metallothionein-III) (MT-III), MT3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MT3