



EPCAM Antibody, Biotin conjugated

Product Code	CSB-PA10799D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16422
Immunogen	Recombinant Human Epithelial cell adhesion molecule protein (24-265AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Epithelial cell adhesion molecule (Ep-CAM) (Adenocarcinoma-associated antigen) (Cell surface glycoprotein Trop-1) (Epithelial cell surface antigen) (Epithelial glycoprotein) (EGP) (Epithelial glycoprotein 314) (EGP314) (hEGP314) (KS 1/4 antigen) (KSA) (Major gastrointestinal tumor-associated protein GA733-2) (Tumor-associated calcium signal transducer 1) (CD antigen CD326), EPCAM, GA733-2 M1S2 M4S1 MIC18 TACSTD1 TROP1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	EPCAM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.