




FTL3 Antibody

Product Code	CSB-PA21589A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	E3NYP3
Immunogen	Recombinant Soybean FT-like protein (1-176AA)
Raised In	Rabbit
Species Reactivity	Soybean
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	FT-like protein, FT2a, FT3 protein, Flowering locus T, Flowering locus T-like protein 3, FTL3, FT, FT2a, FT3, GmFT2a, LOC100814951, GLYMA_16G150700
Immunogen Species	Soybean
Research Area	Others
Target Names	FTL3

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: FTL3 antibody at 5.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 50 kDa

Observed band size: 50 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.