



# FTL3 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA21589C0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	E3NYP3
<b>Immunogen</b>	Recombinant Soybean FT-like protein (1-176AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Soybean
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	FT-like protein, FT2a, FT3 protein, Flowering locus T, Flowering locus T-like protein 3, FTL3, FT, FT2a, FT3, GmFT2a, LOC100814951, GLYMA_16G150700
<b>Immunogen Species</b>	Soybean
<b>Research Area</b>	Others
<b>Target Names</b>	FTL3
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.