



Da-3 Antibody

Product Code	CSB-PA22309A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9SSZ0
Immunogen	Recombinant Trifolium repens DA-3 protein (1-93AA)
Raised In	Rabbit
Species Reactivity	Trifolium repens
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DA-3, Da-3
Immunogen Species	Trifolium repens
Research Area	Others
Target Names	Da-3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.