



LZTS1 Antibody, Biotin conjugated

Product Code	CSB-PA22699D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y250
Immunogen	Recombinant Human Leucine zipper putative tumor suppressor 1 protein (47-212AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine zipper putative tumor suppressor 1 (F37/esophageal cancer-related gene-coding leucine-zipper motif) (Fez1), LZTS1, FEZ1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LZTS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.