



MYH4 Antibody, HRP conjugated

Product Code	CSB-PA22719B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y623
Immunogen	Recombinant Human Myosin-4 protein (896-1045AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin-4 (Myosin heavy chain 2b) (MyHC-2b) (Myosin heavy chain 4) (Myosin heavy chain IIb) (MyHC-IIb) (Myosin heavy chain, skeletal muscle, fetal), MYH4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYH4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.