



OCA2 Antibody, HRP conjugated

Product Code	CSB-PA22759B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q04671
Immunogen	Recombinant Human P protein (2-136AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	P protein (Melanocyte-specific transporter protein) (Pink-eyed dilution protein homolog), OCA2, D15S12 P
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OCA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.