



## NT5C2 Antibody, Biotin conjugated

|                            |                                                                                                      |
|----------------------------|------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA23789D0Rb                                                                                      |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                        |
| <b>Uniprot No.</b>         | P49902                                                                                               |
| <b>Immunogen</b>           | Recombinant Human Cytosolic purine 5'-nucleotidase protein (487-553AA)                               |
| <b>Raised In</b>           | Rabbit                                                                                               |
| <b>Species Reactivity</b>  | Human                                                                                                |
| <b>Tested Applications</b> | ELISA                                                                                                |
| <b>Form</b>                | Liquid                                                                                               |
| <b>Conjugate</b>           | Biotin                                                                                               |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                     |
| <b>Purification Method</b> | >95%, Protein G purified                                                                             |
| <b>Isotype</b>             | IgG                                                                                                  |
| <b>Clonality</b>           | Polyclonal                                                                                           |
| <b>Alias</b>               | Cytosolic purine 5'-nucleotidase (EC 3.1.3.5) (Cytosolic 5'-nucleotidase II), NT5C2, NT5B NT5CP PNT5 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                 |
| <b>Research Area</b>       | Epigenetics and Nuclear Signaling                                                                    |
| <b>Target Names</b>        | NT5C2                                                                                                |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                          |