



gata1-a Antibody

Product Code	CSB-PA326431LA01XBE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23767
Immunogen	Recombinant Xenopus laevis GATA-binding factor 1-A protein (1-359AA)
Raised In	Rabbit
Species Reactivity	Xenopus laevis
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GATA-binding factor 1-A (Transcription factor xGATA-1A), gata1-a, gata1a
Immunogen Species	Xenopus laevis (African clawed frog)
Research Area	Others
Target Names	gata1-a
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.