



TLE4 Antibody, Biotin conjugated

Product Code	CSB-PA523988LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q04727
Immunogen	Recombinant Human Transducin-like enhancer protein 4 protein (284-394AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transducin-like enhancer protein 4 (Grg-4) (Groucho-related protein 4), TLE4, GRG4 KIAA1261
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TLE4