



TLE4 Antibody, Biotin conjugated

Product Code	CSB-PA523988LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q04727
Immunogen	Recombinant Human Transducin-like enhancer protein 4 protein (284-394AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transducin-like enhancer protein 4 (Grg-4) (Groucho-related protein 4), TLE4, GRG4 KIAA1261
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TLE4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.