



PDE9A Antibody, HRP conjugated

Product Code	CSB-PA527985LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76083
Immunogen	Recombinant Human High affinity cGMP-specific 3',5'-cyclic phosphodiesterase 9A protein (426-533AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	High affinity cGMP-specific 3',5'-cyclic phosphodiesterase 9A (EC 3.1.4.35), PDE9A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDE9A