



PDE9A Antibody, Biotin conjugated

Product Code	CSB-PA527985LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76083
Immunogen	Recombinant Human High affinity cGMP-specific 3',5'-cyclic phosphodiesterase 9A protein (426-533AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	High affinity cGMP-specific 3',5'-cyclic phosphodiesterase 9A (EC 3.1.4.35), PDE9A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDE9A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.