



KCNJ8 Antibody, Biotin conjugated

Product Code	CSB-PA613604LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15842
Immunogen	Recombinant Human ATP-sensitive inward rectifier potassium channel 8 protein (299-424AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-sensitive inward rectifier potassium channel 8 (Inward rectifier K(+) channel Kir6.1) (Potassium channel, inwardly rectifying subfamily J member 8) (uKATP-1), KCNJ8
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	KCNJ8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.