



PFKFB4 Antibody, HRP conjugated

Product Code	CSB-PA613704LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16877
Immunogen	Recombinant Human 6-phosphofructo-2-kinase/fructose-2,6-bisphosphatase 4 protein (266-368AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	6-phosphofructo-2-kinase/fructose-2,6-bisphosphatase 4 (6PF-2-K/Fru-2,6-P2ase 4) (PFK/FBPase 4) (6PF-2-K/Fru-2,6-P2ase testis-type isozyme) [Includes: 6-phosphofructo-2-kinase (EC 2.7.1.105); Fructose-2,6-bisphosphatase (EC 3.1.3.46)], PFKFB4
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PFKFB4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.