



PTPN13 Antibody, Biotin conjugated

Product Code	CSB-PA614263LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12923
Immunogen	Recombinant Human Tyrosine-protein phosphatase non-receptor type 13 protein (1965-2173AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein phosphatase non-receptor type 13 (EC 3.1.3.48) (Fas-associated protein-tyrosine phosphatase 1) (FAP-1) (PTP-BAS) (Protein-tyrosine phosphatase 1E) (PTP-E1) (hPTPE1) (Protein-tyrosine phosphatase PTPL1), PTPN13, PNP1 PTP1E PTPL1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PTPN13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.