



PON3 Antibody, HRP conjugated

Product Code	CSB-PA614883LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15166
Immunogen	Recombinant Human Serum paraoxonase/lactonase 3 protein (269-335AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum paraoxonase/lactonase 3 (EC 3.1.1.2) (EC 3.1.1.81) (EC 3.1.8.1), PON3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PON3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.