



SLC14A1 Antibody, Biotin conjugated

Product Code	CSB-PA615678LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13336
Immunogen	Recombinant Human Urea transporter 1 protein (1-51AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Urea transporter 1 (Solute carrier family 14 member 1) (Urea transporter, erythrocyte), SLC14A1, HUT11 JK RACH1 UT1 UTE
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC14A1